



**Orto Specialized**

GC Orthodontics exclusivist  
Global experts in Orthodontics

Engineering By  TOMY

**Orto Specialized**

Carrier de Tarragona 157  
Torre NN - 08014 Barcelona  
info@ortospecialized.es  
www.ortospecialized.es

La técnica  
**BIDIMENSIONALE**  
en una explosión  
de colores nunca  
antes vistos



## LEGEND MEDIUM VERTICAL SLOT

Maxilar	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
U 1	0,41	12°	5°	0°	18x25	32-1110-0034	32-1210-0034
U 2	0,93	8°	9°	0°	18x25	32-1120-0034	32-1220-0034
U 3 hook	0,36	0°	7°	0°	22x28	32-213D-0040	32-223D-0040
U 4/5	0,43	0°	0°	0°	22x28	32-2540-0040	Universal

Mandíbula	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
L 1/2	1,09	0°	0°	0°	18x25	32-1610-0034	Universal
L 3 hook	0,36	0°	5°	0°	22x28	32-243D-0040	32-233D-0040
L 4/5	0,36	0°	0°	0°	22x28	32-2640-0040	Universal

### OPZIONI

U 4/5 hook	0,43	0°	0°	0°	22x28	32-214M-0040	32-224M-0040
L 4/5 hook	0,36	0°	0°	0°	22x28	32-244D-0040	32-234D-0040



## LEGEND MINI VERTICAL SLOT

Maxilar	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
U 1	0,70	12°	5°	0°	18x25	31-1110-0034	31-1210-0034
U 2	0,97	8°	9°	0°	18x25	31-1120-0034	31-1220-0034
U 3 hook	0,56	0°	7°	0°	22x28	31-213D-0040	31-223D-0040
U 4/5	0,56	0°	0°	0°	22x28	31-2540-0040	Universal

Mandíbula	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
L 1/2	0,97	0°	0°	0°	18x25	31-1610-0034	Universal
L 3 hook	0,56	0°	5°	0°	22x28	31-243D-0040	31-233D-0040
L 4/5	0,56	0°	0°	0°	22x28	31-2640-0040	Universal

### OPZIONI

U 4/5 hook	0,56	0°	0°	0°	22x28	31-215M-0040	31-225M-0040
L 4/5 hook	0,56	0°	0°	0°	22x28	31-244D-0040	31-234D-0040



## CHIC

Maxilar	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
U 1	0,5	12°	5°	0°	18x25	14-1110-0010	14-1210-0010
U 2	1	8°	9°	0°	18x25	14-1120-0010	14-1220-0010
U 3 hook	0,42	-2	13°	0°	22x28	14-213D-0010	14-223D-0010
U 4/5	0,6	-7	0°	0°	22x28	14-2540-0010	Universal

Mandíbula	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
L 1/2	0,95	0°	0°	0°	18x25	14-1610-0010	Universal
L 3 hook	0,36	-11°	7°	0°	22x28	14-243D-0010	14-233D-0010



## Self Ligating EXPERIENCE

Maxilar	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
U 1	0,81	12°	5°	0°	18x25	20-1110-0010	20-1210-0010
U 2	1,07	8°	9°	0°	18x25	20-1120-0010	20-1220-0010
U 3 hook	0,46	0°	8°	0°	22x28	20-213D-0030	20-223D-0030
U 4/5	0,46	-7	0°	0°	22x28	20-2140-0010	20-2240-0010

Mandíbula	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
L 1/2	1,17	-1	2°	0°	18x25	20-1410-0010	20-1310-0010
L 3 hook	0,56	0°	3°	0°	22x28	20-243D-0021	20-233D-0021
L 4/5	0,33	-7	-1	0°	22x28	20-2440-0030	20-2340-0030

### OPCIONES

U 4 hook	0,46	-7	0°	0°	22x28	20-214D-0020	20-224D-0020
U 5 hook	0,46	-7	0°	0°	22x28	20-215D-0020	20-225D-0020
L 4/5 hook	0,33	-7	-1	0°	22x28	20-244D-0030	20-234D-0030



## Self Ligating EXPERIENCE MINI

Maxilar	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
U 1	0,74	12°	5°	0°	18x25	21-1110-0010	21-1210-0010
U 2	1,06	8°	9°	0°	18x25	21-1120-0010	21-1220-0010
U 3 hook	0,55	-2	-13	0°	22x28	21-213D-0010	21-223D-0010
U 4/5	0,70	-7	0°	0°	22x28	21-2140-0010	21-2240-0010

Mandíbula	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
L 1/2	1,05	-1	2°	0°	18x25	21-1410-0010	21-1310-0010
L 3 hook	0,55	-6	3°	0°	22x28	21-243D-0020	21-233D-0020
L 4	0,55	-12	2°	0°	22x28	21-2440-0020	21-2340-0020
L 5	0,55	-17	2°	0°	22x28	21-2450-0020	21-2350-0020

### OPCIONES

U 3 passive hook		-2	-13	0°	22x28	21-413D-0010	21-423D-0010
U 4/5 hook		-7	0°	0°	22x28	21-214D-0010	20-224D-0010
U 4/5 passive		-7	0°	0°	22x28	21-4140-0010	21-4240-0010
U 4/5 passive hook		-7	0°	0°	22x28	21-414D-0010	21-424D-0010
L 3 passive hook		-6	3°	0°	22x28	21-443D-0020	21-433D-0020
L 4 passive hook		-12	2°	0°	22x28	21-444D-0020	21-443D-0020
L 4 hook		-12	2°	0°	22x28	21-244D-0020	21-243D-0020
L 4 passive		-12	2°	0°	22x28	21-4440-0020	21-4430-0020
L 5 passive hook		-17	2°	0°	22x28	21-445D-0020	21-435D-0020
L 5 hook		-17	2°	0°	22x28	21-245D-0020	21-235D-0020
L 5 passive		-17	2°	0°	22x28	21-4450-0020	21-4350-0020



## Self Ligating EXPERIENCE MINI Rhodium

Maxilar	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
U 1	0,74	12°	5°	0°	18x25	24-1110-0010	24-1210-0010
U 2	1,06	8°	9°	0°	18x25	24-1120-0010	24-1220-0010
U 3 hook	0,55	-2	-13	0°	22x28	24-213D-0010	24-223D-0010
U 4/5	0,70	-7	0°	0°	22x28	24-2140-0010	24-2240-0010

Mandíbula	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
L 1/2	1,05	-1	2°	0°	18x25	24-1410-0010	24-1310-0010
L 3 hook	0,55	-11	7°	0°	22x28	24-243D-0010	24-233D-0010
L 4	0,55	-17	1	0°	22x28	24-2440-0010	24-2340-0010
L 5	0,55	-22	1	0°	22x28	24-2450-0010	24-2350-0010

### OPCIONES

U 4/5 hook		-7	0°	0°	22x28	24-214D-0010	24-224D-0010
L 4 hook		-17	1	0°	22x28	24-244D-0010	24-234D-0010
L 5 hook		-22	1	0°	22x28	24-245D-0010	24-235D-0010



## Self Ligating EXPERIENCE CLEAR

Maxilar	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
U 1	0,81	12°	5°	0°	18x25	10-1110-0010	10-1210-0010
U 2	1,07	8°	9°	0°	18x25	10-1120-0010	10-1220-0010
U 3 hook	0,46	0°	8°	2°	22x28	10-213U-0030	10-223U-0030
U 4/5	0,46	-7	0°	2°	22x28	10-2140-0010	10-2240-0010

Mandíbula	In out	Torque	Ang.	Offset	Slot size	Derecha	Izquierda
L 1/2	1,17	1°	0°	0°	18x25	10-1610-0030	Universale
L 3 hook	0,56	0°	3°	0°	22x28	10-243U-0021	10-233U-0021
L 4	0,38	-17	0°	4°	22x28	10-2450-Z010	10-2350-Z010
L 5	0,28	-22	0°	4°	22x28	10-2450-Z010	10-2350-Z010

### OPCIONES

U 4/5 hook	0,46	-7	0°	2°	22x28	10-214U-0010	10-224U-0010
L 4 hook	0,38	-17	0°	4°	22x28	10-244U-Z010	10-234U-Z010
L 5 hook	0,28	-22	0°	4°	22x28	10-245U-Z010	10-235U-Z010

With Rhodium coating



Without Rhodium coating



## B-Sure by AMG

### Experience Mini Rhodium



## Initialloy round archwires

Fuerza	Dimensiones	Large P/N		Medium P/N	
		Maxilar	Mandíbula	Maxilar	Mandíbula
Light	.014"	71-A1U0-0014	71-A1L0-0014	-	-
Light	.016"	71-A1U0-0016	71-A1L0-0016	-	-
Light	.018"	71-A1U0-0018	71-A1L0-0018	-	-
Medium	.012"	71-A2U0-0012	71-A2L0-0012	-	-
Medium	.014"	71-A2U0-0014	71-A2L0-0014	71-B2U0-0014	71-B2L0-0014
Medium	.014" (Dimple)	71-A2UD-0014	71-A2LD-0014	71-B2UD-0014	71-B2LD-0014
Medium	.016"	71-A2U0-0016	71-A2L0-0016	71-B2U0-0016	71-B2L0-0016
Medium	.016" (Dimple)	71-A2UD-0016	71-A2LD-0016	71-B2UD-0016	71-B2LD-0016
Medium	.018"	71-A2U0-0018	71-A2L0-0018	71-B2U0-0018	71-B2L0-0018
Medium	.018" (Dimple)	71-A2UD-0018	71-A2LD-0018	71-B2UD-0018	71-B2LD-0018
Medium	.020"	71-A2U0-0020	71-A2L0-0020	-	-
Heavy	.014"	71-A3U0-0014	71-A3L0-0014	-	-
Heavy	.016"	71-A3U0-0016	71-A3L0-0016	-	-
Heavy	.018"	71-A3U0-0018	71-A3L0-0018	-	-

Fuerza	Dimensio	Universal P/N	
		Maxilar	Mandíbula
Light	.014"	71-C1U0-0014	71-C1L0-0014
Light	.016"	71-C1U0-0016	71-C1L0-0016
Light	.018"	71-C1U0-0018	71-C1L0-0018
Medium	.014"	71-C2U0-0014	71-C2L0-0014
Medium	.014" (Dimple)	71-C2UD-0014	71-C2LD-0014
Medium	.016"	71-C2U0-0016	71-C2L0-0016
Medium	.016" (Dimple)	71-C2UD-0016	71-C2LD-0016
Medium	.018"	71-C2U0-0018	71-C2L0-0018
Medium	.018" (Dimple)	71-C2UD-0018	71-C2LD-0018
Heavy	.016"	71-C3U0-0016	71-C3L0-0016



## BioEdge

Fuerza	Dimensiones	Large P/N		Medium P/N	
		Maxilar	Mandíbula	Maxilar	Mandíbula
<b>80g.</b>					
	.016" x .016"	72-A1U0-1616	72-A1L0-1616	-	-
	.016" x .022"	72-A1U0-1622	72-A1L0-1622	-	-
<b>100g.</b>					
	.020" x .020"	72-A1U0-2020	72-A1L0-2020	-	-

Fuerza Dimensiones	Large P/N		Medium P/N	
	Maxilar	Mandíbula	Maxilar	Mandíbula
<b>160g.</b>				
.016"x.016"	72-A2U0-1616	72-A2L0-1616	-	-
.016"x.016" (Dimple)	72-A2UD-1616	72-A2LD-1616	72-B2UD-1616	72-B2LD-1616
.016"x.022"	72-A2U0-1622	72-A2L0-1622	72-B2U0-1622	72-B2LO-1622
.016"x.022" (Dimple)	72-A2UD-1622	72-A2LD-1622	72-B2UD-1622	72-B2LD-1622
<b>200g.</b>				
.017"x.025"	72-A2U0-1725	72-A2L0-1725	-	-
.017"x.025" (Dimple)	72-A2UD-1725	72-A2LD-1725	72-B2UD-1725	72-B2LD-1725
.018"x.018"	72-A2U0-1818	72-A2L0-1818	-	-
.018"x.025"	72-A2U0-1825	72-A2L0-1825	72-B2U0-1825	72-B2LO-1825
.019"x.025"	72-A2U0-1925	72-A2L0-1925	-	-
.020"x.020"	72-A2U0-2020	72-A2L0-2020	-	-
.021"x.028"	72-A2U0-2128	72-A2L0-2128	-	-
<b>240g.</b>				
.016"x.022"	72-A3U0-1622	72-A3L0-1622	-	-
<b>300g.</b>				
.018"x.025"	72-A3U0-1825	72-A3L0-1825	-	-

Fuerza Dimensiones	Ovoid (Brader) P/N	
	Maxilar	Mandíbula
<b>80g.</b>		
.016"x.016"	72-C1U0-1616	72-C1L0-1616
.016"x.022"	72-C1U0-1622	72-C1L0-1622
<b>100g.</b>		
.017"x.025"	72-C1U0-1725	72-C1L0-1725
.018"x.018"	72-C1U0-1818	72-C1L0-1818
.018"x.025"	72-C1U0-1825	72-C1L0-1825
.020"x.020"	72-C1U0-2020	72-C1L0-2020
<b>160g.</b>		
.016"x.016"	72-C2U0-1616	72-C2L0-1616
.016"x.016" (Dimple)	72-C2UD-1616	72-C2LD-1616
.016"x.022"	72-C2U0-1622	72-C2L0-1622
.016"x.022" (Dimple)	72-C2UD-1622	72-C2LD-1622
<b>200g.</b>		
.017"x.025"	72-C2U0-1725	72-C2L0-1725
.017"x.025" (Dimple)	72-C2UD-1725	72-C2LD-1725
.018"x.018"	72-C2U0-1818	72-C2L0-1818
.018"x.025"	72-C2U0-1825	72-C2L0-1825
.020"x.020"	72-C2U0-2020	72-C2L0-2020
<b>200g.</b>		
.016"x.022"	72-C3U0-1722	72-C3L0-1722

Fuerza Dimensiones	Standart	
	Maxilar	Mandíbula
<b>80g.</b>		
.016"x.022"	72-E1U0-1622	72-E1L0-1622
<b>100g.</b>		
.018"x.025"	72-E1U0-1825	72-E1L0-1825

Fuerza Dimensiones	Universali P/N	
	Maxilar	Mandíbula
<b>160g.</b>		
.014"x.025"	72-D200-1425	
.016"x.022"	72-D200-1622	
<b>200g.</b>		
.018"x.025"	72-D200-1825	
.020"x.020"	72-D200-2020	



## BioActive

Dimensiones	Large P/N		Medium P/N	
	Maxilar	Mandíbula	Maxilar	Mandíbula
.016"x.016"	70-A0U0-1616	70-A0L0-1616	-	-
.016"x.022"	70-A0U0-1622	70-A0L0-1622	-	-
.016"x.022" (Dimple)	70-A0UD-1622	70-A0LD-1622	70-B0U0-1622	70-B0L0-1622
.018"x.018"	70-A0U0-1818	70-A0L0-1818	70-B0U0-1818	70-B0L0-1818
.018"x.018" (Dimple)	70-A0UD-1818	70-A0LD-1818	70-B0UD-1818	70-B0LD-1818
.018"x.025"	70-A0U0-1825	70-A0L0-1825	70-B0U0-1825	70-B0L0-1825
.019"x.025"	70-A0U0-1925	70-A0L0-1925	70-B0U0-1925	70-B0L0-1925
.020"x.020"	70-A0U0-2020	70-A0L0-2020	70-B0U0-2020	70-B0L0-2020
.020"x.020" (Dimple)	70-A0UD-2020	70-A0LD-2020	70-B0UD-2020	70-B0LD-2020

Dimensiones	Ovoid (Brader) P/N	
	Maxilar	Mandíbula
.016"x.016"	70-C0U0-1616	70-C0L0-1616
.016"x.022"	70-C0U0-1622	70-C0L0-1622
.016"x.022" (Dimple)	70-C0UD-1622	70-C0LD-1622
.018"x.018"	70-C0U0-1818	70-C0L0-1818
.018"x.018" (Dimple)	70-C0UD-1818	70-C0LD-1818
.018"x.025"	70-C0U0-1825	70-C0L0-1825
.019"x.025"	70-C0U0-1925	70-C0L0-1925
.020"x.020"	70-C0U0-2020	70-C0L0-2020
.020"x.020" (Dimple)	70-C0UD-2020	70-C0LD-2020



## Molybdenum Beta Titanium

Dimensiones	Large P/N		Medium P/N	
	Maxilar	Mandíbula	Maxilar	Mandíbula
.016"	73-A0U0-0016	73-A0L0-0016	-	-
.018"	73-A0U0-0018	73-A0L0-0018	-	-
.016"x.022"	73-A0U0-1622	73-A0L0-1622	73-B0U0-1622	73-B0L0-1622
.017"x.025"	73-A0U0-1725	73-A0L0-1725	73-B0U0-1725	73-B0L0-1725
.018"x.018"	73-A0U0-1818	73-A0L0-1818	73-B0U0-1818	73-B0L0-1818
.018"x.025"	73-A0U0-1825	73-A0L0-1825	73-B0U0-1825	73-B0L0-1825
.019"x.025"	73-A0U0-1925	73-A0L0-1925	73-B0U0-1925	73-B0L0-1925
.021"x.025"	73-A0U0-2125	73-A0L0-2125	-	-

Dimensiones	Ovoid (Brader) P/N	
	Maxilar	Mandíbula
.016"x.022"	73-C0U0-1622	73-C0L0-1622
.017"x.025"	73-C0U0-1725	73-C0L0-1725
.018"x.018"	73-C0U0-1818	73-C0L0-1818
.018"x.025"	73-C0U0-1825	73-C0L0-1825
.019"x.025"	73-C0U0-1925	73-C0L0-1925
.021"x.025"	73-C0U0-2125	73-C0L0-2125

## Monostranded stainless steel

Dimensiones	Large P/N ST		Medium P/N ST	
	Maxilar	Mandíbula	Maxilar	Mandíbula
.014"	78-A2U0-0014	78-A2L0-0014	-	-
.016"	78-A2U0-0016	78-A2L0-0016	78-B2U0-0016	78-B2L0-0016
.018"	78-A2U0-0018	78-A2L0-0018	78-B2U0-0018	78-B2L0-0018
.020"	78-A2U0-0020	78-A2L0-0020	-	-
.016"x.016"	78-A2U0-1616	78-A2L0-1616	-	-
.016"x.022"	78-A2U0-1622	78-A2L0-1622	78-B2U0-1622	78-B2L0-1622
.017"x.025"	78-A2U0-1725	78-A2L0-1725	78-B2U0-1725	78-B2L0-1725
.018"x.018"	78-A2U0-1818	78-A2L0-1818	78-B2U0-1818	78-B2L0-1818
.018"x.022"	78-A2U0-1822	78-A2L0-1822	-	-
.018"x.025"	78-A2U0-1825	78-A2L0-1825	78-B2U0-1825	78-B2L0-1825
.019"x.019"	78-A2U0-1919	78-A2L0-1919	-	-
.019"x.025"	78-A2U0-1925	78-A2L0-1925	78-B2U0-1925	78-B2L0-1925
.020"x.020"	78-A2U0-2020	78-A2L0-2020	78-B2U0-2020	78-B2L0-2020
.021"x.025"	78-A2U0-2125	78-A2L0-2125	78-B2U0-2125	78-B2L0-2125
.021"x.028"	78-A2U0-2128	78-A2L0-2128	78-B2U0-2128	78-B2L0-2128

Dimensiones	Large P/N ST		Medium P/N SHT	
	Maxilar	Mandíbula	Maxilar	Mandíbula
.014"	78-C2U0-0014	78-C2L0-0014	-	-
.016"	78-C2U0-0016	78-C2L0-0016	-	-
.018"	78-C2U0-0018	78-C2L0-0018	-	-
.020"	78-C2U0-0020	78-C2L0-0020	-	-
.016"x.016"	78-C2U0-1616	78-C2L0-1616	78-C3U0-1616	78-C3L0-1616
.016"x.022"	78-C2U0-1622	78-C2L0-1622	78-C3U0-1622	78-C3L0-1622
.017"x.017"	78-C2U0-1717	78-C2L0-1717	-	-
.017"x.025"	78-C2U0-1725	78-C2L0-1725	-	-
.018"x.022"	-	-	78-C3U0-1825	78-C3L0-1825
.018"x.025"	78-C2U0-1825	78-C2L0-1825	-	-
.019"x.025"	78-C2U0-1925	78-C2L0-1925	-	-
.020"x.020"	78-C2U0-2020	78-C2L0-2020	-	-
.020"x.025"	78-C2U0-2025	78-C2L0-2025	-	-
.021"x.025"	78-C2U0-2125	78-C2L0-2125	-	-

Dimensiones	Standart	
	Maxilar	Mandíbula
<b>Super Temperato Gold</b>		
.016"x.022"	78-E5U0-1622	78-E5L0-1622
<b>Super Temperato Gold</b>		
.018"x.022"	78-E5U0-1822	78-E5L0-1822





## Banda con tubos doble no convertible Bidimensionale

I° molar superior	Codice
Tubo Doble gingival o Oclusal	41-2160/2260-0080 + Cleat Linguale 43-00220-0000
I° molar inferior	Codice
Tubo Doble gingival	41-2160/2260-0080 + Cleat Linguale 43-00220-0000

## Banda con tubos doble convertible Bidimensionale

I° Molare superiore	Codice
Tubo Doble gingival o Oclusal	41-2160/2260-0082 + Cleat Linguale 43-00220-0000
I° molar inferior	Codice
Tubo Doble gingival	41-2160/2260-0082 + Cleat Linguale 43-00220-0000



## Banda con tubos con una sola canula Bidimensionale

II° molar superior	Codice
Tubo con una sola cánula	41- 2370/2470- 0096 + cleat linguale 43-0220-0000
II° molar inferior	Codice
Tubo con una sola cánula	41- 2370/2470- 0096 + cleat linguale 43-0220-0000